

Omark Life Sciences provides RODAC (Replicate Qrganism Detection and Counting) plates for monitoring surface contamination for various purposes. Plates are made with high quality polystyrene

Section 1 : Gridded and Convex bottom 65*12mm

Section 2 : Gamma Sterile Disposable Scored Top Petri Plates

Section 3: Hazards Identification Hazard: Not classified as hazardous.

Section 4: First Aid Measures: No specific measures necessary.

Section 5: Fire Fighting Measures: Not combustible.

Section 6 : Accidental Release Measures : No specific measures necessary.

Section 7 : Storage - Store below 30°C

Section 8: Exposure Controls / Personal Protection

Wear appropriate NIOSH/MSHA-approved respirator, chemical-resistant gloves, safety goggles, other protective clothing. Mechanical exhaust required.

Section 9 : Physical and Chemical Properties

Appearance: Sterile Disposable petri plates, optically clear, 65 mm diameter X 12 mm

Section 10: Stability and Reactivity Stability:

Product is stable if stored as per the conditions specified under storage below 30°C

Product loses its potency/performance above 45°C. Conditions to avoid: Heat and light.

Section 11: Toxicological Information Non toxic.

